



**FARETEC INC**  
RESPONDER PRODUCTS



**CT-7 LEG TRACTION SPLINT INSTRUCTIONS  
FOR ADULT, PEDIATRIC AND BELOW KNEE APPLICATIONS**

1.)



Open zipper. Remove ankle hitch from bag, this will be applied first.

Remove splint from bag. Hold upright and allow tubes to intersect. Lay alongside patients injured leg.

Top of splint should be near the patient's hip crest.

Remove Ischial Cap and disengage tubing as necessary in order to shorten the splint. Fold excess tubes alongside and replace cap.

2.)



3.)



4.)

On the ankle hitch, make sure the buckle is pointed toward patient's foot.

Wrap **TIGHTLY** around base of ankle. Buckle should be on **OUTSIDE** of leg.

5.)



Attached buckle at bottom of splint to the ankle hitch.

Unclip one side of the Ischial Strap and wrap around the top of the patients inner thigh. Clip buckle back onto splint.

6.)





**FARETEC INC**  
RESPONDER PRODUCTS



**CT-7 LEG TRACTION SPLINT INSTRUCTIONS  
FOR ADULT, PEDIATRIC AND BELOW KNEE APPLICATIONS**



7.)

Locate tail of Ischial Strap and pull **TIGHT**. It should fit snug around patients leg.



8.)

Attach the leg straps.

Wrap one **BELOW** the knee.

Wrap the second leg strap **ABOVE** the knee but not directly over the injury.



9.)

(Extension strap available for wrap around both patients legs for increased stability & immobilization.)



10.)

Check distal pulse and other possible contradictions for leg traction.

Spin traction mechanism until appropriate traction is applied.

- Legs should be equal length.
- Reduce swelling at point of injury
- Increase patient comfort/relief.

Works equally with below the knee application . Remove ischial cap and shorten the splint accordingly. Wrap ankle hitch tightly above the patients calf.

